

PharmaOne AI – UI/UX Requirements (Updated)

Authentication, demo-user admin flow, Deviation intake, AI review and Reports.

1. Scope

PharmaOne AI is an AI-powered Deviation Intake and Management workflow for an API pharmaceutical manufacturer. Core flow: deviation document/text → AI extraction → Log Deviation form → AI impact/severity assistance → user review/edit → save → reports.

Customer Complaint and other QMS modules are not implemented as full modules. Company switching is intentionally excluded from this version.

2. Authentication and Access

Use one login page for Admin and normal users with Supabase Auth. No public Sign Up or self-service company registration. One Admin is manually provisioned. Admin creates demo users and assigns company, site and one of three roles: QA Manager, Production User, QC User.

Each user has one company membership. No company selector or company switching.

3. Navigation

Navigation	Purpose
Dashboard	Overview of deviation activity and key metrics.
Deviations	Create, process, review and view deviations.
Reports	View saved deviation reporting.
Notifications	Notifications area.
Profile	User profile and logout.
Admin Panel	Admin-only demo user management and audit logs.

4. Deviation Intake UX

User uploads a deviation PDF or provides text/email content. AI processes it and populates the Log Deviation form.

Show processing states: Uploading → Extracting/OCR if required → AI extraction → Knowledge retrieval → Risk assessment → Ready for review.

Do not require manual completion of the whole form before AI processing.

5. Form Fields

Group	Fields
Identification	Deviation ID, date/time, reported by, source
Organization	Company, site, department, responsible team
Product	Product, product code, batch
Event	Manufacturing stage, equipment, title, description
Conditions	Expected condition, actual condition, duration, parameter
Immediate response	Immediate action, batch status, QA/supervisor notification
Assessment	Impact, severity, reason

6. AI Review

Display extracted information and AI recommended impact/severity with a short reason and supporting evidence. User can edit AI-populated fields before saving. Final saved values are the user's reviewed values.

7. Admin UX

Admin Panel contains Demo Users and Audit Logs. Create Demo User fields: full name, work email, employee ID, department, job title, company, site, role and status. Admin cannot create another Admin from this form.

8. Errors

Handle invalid login, inactive user, missing membership, invalid PDF, OCR failure, AI failure, invalid structured output, database failure and session expiry. Never save fabricated AI output.